

C.A.R.

PRODUCTS[®]

makers of fine cleaners

630 Beaulieu Street, Holyoke, MA 01040
Phone 1-800-537-7797 Fax 1-413-536-9979

**GRAPE SEALER
PART #ICS-716
TITRATION PROCEDURE
11/1/09
DRYING AGENT/SEALANT/PROTECTANT KIT**

1. RINSE AND FILL 25 ML SAMPLE TUBE (#9198) TO THE 10 ML MARK WITH WATER TO BE TESTED.
2. USING 1.0 ML PIPET (#4030), ADD 0.5 ML COMPLEXING REAGENT (R-0950). SWIRL TO MIX.
3. ADD 2 DROPS R-0881 TOLUIDINE BLUE O INDICATOR. SWIRL TO MIX. SAMPLE SHOULD BE DARK BLUE.
4. ADD R-0884 QAC TITRATING SOLUTION DROPWISE, SWIRLING AND COUNTING AFTER EACH DROP, UNTIL COLOR CHANGES FROM DARK BLUE TO PURPLE. ALWAYS HOLD BOTTLE IN VERTICAL POSITION.

<u>USAGE</u>	<u>DILUTION</u>	<u>DROPS</u>
½ FL. OZ.	1:256	2-3
*1 FL. OZ.	1:128	5
1.5 FL. OZ.	1:96	7.8
2 FL. OZ.	1:64	10
* RECOMMENDED RANGE		